

**THE TECHNOLOGY YOU WANT.
THE FLEXIBILITY YOU NEED.**

The advanced technology of Vbeam is more accessible than ever thanks to Pathways—your path to unprecedented Choice, Flexibility and Value. The **choice** of a complete range of innovative technologies for optimal treatment.

The **flexibility** to meet your specific needs today, then easily upgrade as you expand your practice.

The **value** that comes from cost-effectively growing your business—step-by-step, with the leader in medical and aesthetic lasers.

To learn more about Pathways, visit Candela at www.candelalaser.com/pathways.cfm or call 1-800-821-2013.



CHOICE FLEXIBILITY VALUE



CE 0123

Complies with FDA performance standards for laser products except for deviations pursuant to laser Notice No. 50 dated July 26, 2002

**Patent pending.

Specifications subject to change without notice.

Vbeam is a registered trademark of Candela Corporation. Dynamic Cooling Device and DCD are trademarks.

©2005 Candela Corporation. All rights reserved. Printed in the USA. 11/05 0920-20-0053 Rev. 04



Candela Corporation
530 Boston Post Road
Wayland, MA 01778 USA
Phone: (508) 358-7637
Fax: (508) 358-5569
Toll Free: (800) 821-2013
www.candelalaser.com



DCD

THE MOST ADVANCED PULSED DYE LASER TECHNOLOGY



Choose the One That's Right for You



THE MOST TRUSTED PULSED DYE LASER IS ALSO THE MOST ADVANCED

Candela's Vbeam® pulsed dye laser has a 20-year proven track record of safety, reliability and efficacy for treatment of a wide range of vascular conditions. So how did we improve on the gold standard? Three ways: advanced technology; increased versatility; and simplified ease of use.

At Candela, our pursuit of excellence means we never stop refining our products and services. The end result is effective treatments for your patients, backed by Candela's unsurpassed training, marketing and service expertise.


Next Generation Technology, Unprecedented Product Versatility

Candela's gold standard Vbeam laser is better than ever thanks to cutting edge technology now available in three different configurations. One of them is right for your practice. Best of all, like all Candela products, our lasers are designed to grow with your practice. They are fully upgradeable, with training and service designed to keep pace with your changing needs.

A DYNAMIC COOLING DEVICE SO ORIGINAL, IT'S PATENTED

Candela's patented Dynamic Cooling Device™ (DCD™) provides consistent, reproducible epidermal cooling and unparalleled comfort for your patients. DCD prepares skin for treatment by selectively cooling the upper layers with a burst of cryogen milliseconds before the laser pulse. DCD's unique design provides the user with excellent visibility throughout the entire process. The most powerful benefits to you? Speed and Patient Safety. DCD is operator technique independent—it will provide the same epidermal protection no matter how fast the laser is fired. The more patients you see and treat optimally, the more profitable for your practice..

with
DCD



SUPPORT AFTER THE SALE

Any laser purchase includes access to MyCandela.com, an exclusive marketing resource for your new laser. Services include FREE Website Builder, Laser Marketing Tool Kits, Learning center and additional clinical information on the capabilities of your Candela laser.*

*Availability varies by country.

OUR PATENTED PIGMENTED LESION
HANDPIECE** TARGETS AND ELIMINATES
BROWN SPOTS AT THE EPIDERMIS LEVEL
WITHOUT IMPACTING UNDERLYING
BLOOD VESSELS CAUSING PURPURA.

